3M™ Cavityon™ Advanced Hudbeskyttelse

3M Norge

Annelise Vågan Pettersen

3M Sårprodukter - unike løsninger

Unik klebeteknologi

Hypoallergisk

Spesialdesign

Problemløsere

Cavilon™ No Sting barrirefilm

langvarig hudbeskyttelse

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


*Sørger for optimal
pasientpleie*

IAD er et vanlig problem

¹ Bliss et al. J Wound Ostomy Continence Nurs. 2011; ²Kottner et al. Int J Nurs Stud. 2014; ³Campbell et al. Int Wound J. 2014; ⁴Long et al. J Wound Ostomy Continence Nurs. 2012; ⁵Gray M and Bartos S. Presented at the 23rd Annual Meeting of the Wound Healing Society. 2013; ⁶Gray M et al. Journal of Wound, Ostomy and Continence Nursing. 2007; ⁷Nix and Haugen. Drugs Aging 2010;27(6):491-6

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Diagnosis: categorization

Category 1: Persistent redness

1A - Persistent redness without clinical signs of infection

Critical criterion

- Persistent redness
A variety of tones of redness may be present. Patients with darker skin tones, the skin may be paler or darker than normal, or purple in colour.

Additional criteria

- Marked areas or discolouration from a previous (healed) skin defect
- Shiny appearance of the skin
- Macerated skin
- Intact vesicles and bullae
- Skin may feel tense or swollen at palpation
- Burning, tingling, itching or pain

1A

Category 2: Skin loss

2A - Skin loss without clinical signs of infection

Critical criterion

- Skin loss
Skin loss may present as skin erosion (may result from damaged/eroded vesicles or bullae), denudation or excoriation. The skin damage pattern may be diffuse.

Additional criteria

- Persistent redness
A variety of tones of redness may be present. Patients with darker skin tones, the skin may be paler or darker than normal, or purple in colour
- Marked areas or discolouration from a previous (healed) skin defect
- Shiny appearance of the skin
- Macerated skin
- Intact vesicles and bullae
- Skin may feel tense or swollen at palpation
- Burning, tingling, itching or pain

2A

1B - Persistent redness with clinical signs of infection

Critical criteria

- Persistent redness
A variety of tones of redness may be present. Patients with darker skin tones, the skin may be paler or darker than normal, or purple in colour.
- Signs of infection
Such as white scaling of the skin (suggesting a fungal infection) or satellite lesions (pustules surrounding the lesion, suggesting a *Candida albicans* fungal infection).

Additional criteria

- Marked areas or discolouration from a previous (healed) skin defect
- Shiny appearance of the skin
- Macerated skin
- Intact vesicles and bullae
- The skin may feel tense or swollen at palpation
- Burning, tingling, itching or pain

1B

2B - Skin loss with clinical signs of infection

Critical criteria

- Skin loss
Skin loss may present as skin erosion (may result from damaged/eroded vesicles or bullae), denudation or excoriation. The skin damage pattern may be diffuse.
- Signs of infection
Such as white scaling of the skin (suggesting a fungal infection) or satellite lesions (pustules surrounding the lesion, suggesting a *Candida albicans* fungal infection), slough visible in the wound bed (yellow/brown/greyish), green appearance within the wound bed (suggesting a bacterial infection with *Pseudomonas aeruginosa*), excessive exudate levels, purulent exudate (pus) or a shiny appearance of the wound bed.

Additional criteria

- Persistent redness
A variety of tones of redness may be present. Patients with darker skin tones, the skin may be paler or darker than normal, or purple in colour
- Marked areas or discolouration from a previous (healed) skin defect
- Shiny appearance of the skin
- Macerated skin
- Intact vesicles and bullae
- Skin may feel tense or swollen at palpation
- Burning, tingling, itching or pain

2B

CATEGORY 1 = PERSISTENT REDNESS

- CAT 1 A = Persistent redness without clinical signs of infection
- CAT 1 B = Persistent redness with clinical signs of infection

CATEGORY 2 = SKIN LOSS

- CAT 2 A = Skin loss without clinical signs of infection
- CAT 2 B = Skin loss with clinical signs of infection

University Centre for
Nursing and Midwifery

Beeckman D, Van den Bussche K, Beele H, Ciprandi G, Coyer F, De Meyer D, Dunk AM, Fouri A, Garcia-Molina P, Gray M, Johansen E, Karadağ A, Kis Dadara Z, Meaume S, Pokorna A, Romanelli M, Ruppert S, Schoonhoven L, Smet S, Steininger A, Stockmayr M, Van Damme N, Voegeli D, Verhaeghe S, Van Hecke A, Kottner J. (2017). Development and Clinical Validation of the Ghent Global IAD Categorisation tool (GLOBIAD). In progress.

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Medisin

Pasientens kjennetegn

Source: ¹ Beekman et al, Wounds international 2015

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www.3M.no/Cavilon

J428757. 3M and Cavilon are registered trademarks of the 3M company. © 3M 2015. All rights reserved. BR63D88-INC180629

3M

3M

A Case Study using a novel elastomeric barrier for the treatment of severe Incontinence Associated Dermatitis (IAD) in an acute and community setting

Fiona Russell, Nurse Consultant Tissue Viability. NHS Grampian

12/12/2017: Presentation to Tissue Viability Service and 1st application of Cavilon Advanced (in hospital)

15/12/2017: 2nd application of Cavilon Advanced (in hospital)

18/12/2017: 3rd application of Cavilon Advanced (discharged home)

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Patient:

Female patient with IAD for several years, placed in Frederikssund Kommune Denmark. Patient have tested different ointments and barrier creams without any change in results. In november 2017 the patient was introduced to Cavilon™ Advanced Skin Protectant®.

Picture 1: 2017-11-07; picture taken before first application. Treatment start with applications two times per week.

Picture 2: 2017-11-20; ongoing treatment with applications two times per week.

Picture 3: 2017-11-27; ongoing treatment with applications two times per week.

Picture 4: 2018-01-04; patient healed but Cavilon™ Advanced Skin Protectant® is used prophylactic when/if upcoming redness appears again.

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
References

1. O'Meara S, Cullum N, Nelson EA, Dumville JC. Compression for venous leg ulcers. Cochrane Database Syst Rev 2012; 11: CD000265.
2. Nelzen O. Leg ulcers: economic aspects. Phlebology 2000; 15: 110-14.
3. Nelson EA, Bell-Syer SE. Compression for preventing recurrence of venous ulcers. Cochrane Database Syst Rev 2014; 9: CD002303.
4. Sauer K, Rothgang H, Glaeske G. BARMER GEK Heil- und Hilfsmittelreport 2014.
5. Haute Autorite de Sante. Managing venous leg ulcers (excluding dressings). June 2006
6. Deutsche Gesellschaft für Phlebologie (DGP). Guidelines for diagnosis and therapy of venous ulcers (version 8 2008). Phlebologie 2008; 6: 308-29.
7. Gallenkemper G, Wilm S. Leitlinie zu Diagnostik und Therapie des Ulcus cruris venosum der DGP – Kurzversion Herbst 2008. Phlebologie 2010; 5: 290-92.
8. Scottish Intercollegiate Guidelines Network (SIGN). Management of chronic venous leg ulcers. A national clinical guideline. DIGN 2010. European Dermatology Forum (EDF). Guideline on venous leg ulcer version 4.0.
9. European Dermatology Forum (EDF). Guideline on venous leg ulcer version 4.0
10. CONUEI. Conferencia nacional de consenso sobre ulceras de la extremidad inferior: documento de consenso. Barcelona: Edikamed; 2009.
11. Haute Autorite de Sante. Evaluation des dispositifs de compression medicale a usage individuel — Utilisation en pathologies vasculaires. September 2010.
12. NHG Guideline - Venous Ulcers (summary). Dutch College of General Practitioners (NHG), 2010.
13. National Institute for Health and Care Excellence (NICE). Clinical Knowledge Summaries: Leg ulcer - venous 2012.
14. Initiative Chronische Wunden (ICW). Recommendations for compression therapy for patients with venous ulcers. Consensus recommendation. EWMA Journal 2013; 13(2): 42-47.
15. O'Donnell TF, Passman MA, Marston WA, et al. Management of venous leg ulcers: Clinical practice guidelines of the Society for Vascular Surgery and the American Venous Forum. J Vasc Surg 2014; 60: 3S-59S.
16. Petherick ES, Cullum NA, Pickett KE. Investigation of the Effect of Deprivation on the Burden and Management of Venous Leg Ulcers: A Cohort Study Using the THIN Database. PLoS One 2013; 8(3): e58948.
17. Begarin L, Beaujour A, Fainsilber P, et al. [Compression and venous leg ulcer: Observational study in general medicine]. J Mal Vasc 2014; 39(6): 382-88
18. Mosti G, Partsch H. Compression stockings with a negative pressure gradient have a more pronounced effect on venous pumping function than graduated elastic compression stockings. Eur J Vasc Endovasc Surg 2011; 42(2): 261-66.
19. Mosti G, Partsch H. High compression over the calf is more effective than graduated compression in enhancing venous pump function. Eur J Vasc Endovasc Surg 2012; 44: 332-36.